COVID-19 ANTIGEN RAPID TEST (Latex)

JOINSTAR BIOMEDICAL TECHNOLOGY CO., LTD.

CONTACT:

Mobile: +86-13957167265 (Stanley) E-mail: info@joinstar.cn Tel: +86-571-89023160 Fax: +86-571-89028228

Web: <http://en.joinstar.cn>

 @JoinstarBioTech

 @JoinstarBioTech

 www.linkedin.com/company/joinstar

COMPANY PROFILE

Founded in December 2010, with a total investment of 43 million USD, a total construction area of 45929 square meters and a registered capital of 12 million USD, **Joinstar Biomedical Technology Co., Ltd.** is a biomedical technology enterprise focusing on medical in vitro diagnostic products. Joinstar has established a strategic partnership with Morgan Stanley, Including Zhejiang provincial postdoctoral workstation, liquid-phase chip Institute, and Zhejiang provincial centralized supervision platform for Importing & Exporting of special biological products. The company now has a complete set of R & D, production equipment and purification workshop for in vitro diagnostic instruments / reagents / raw materials for POCT, biochemistry, immunity and molecular diagnosis, with an annual production capacity of 20 million diagnostic kits per shift, a 8668 square meters IVD reagent production base and a 1700 square meters public experimental platform. There are 435 employees, accounting for 60% of bachelor degree or above, including nearly 80 doctors, masters and senior titles. 73 patents have been applied for or authorized, including 30 invention patents and 12 international patents. Has obtained 45 NMPA medical device product registration certificates, 29 CE certificates, and more than 40 products under research and registration. The products are sold in more than 50 countries amount the world and 30 provinces, municipalities and autonomous regions in China. It has established a relatively perfect marketing channel with nearly 7000 customers.


JOINSTAR Campus

Workshop


International Cooperation

JOINSTAR & Boditech

JOINSTAR & Axis Shield

JOINSTAR & Skyla

JOINSTAR & Mexican partner

JOINSTAR & Indian partner

JOINSTAR & Omega Diagnostics

Domestic cooperation

JOINSTAR@Wuhan Leishenshan Hospital

Shanghai Jiaotong University-JOINSTAR Joint Laboratory Unveiling

With Professor Sun Baoqing, Guangzhou Institute of Respiratory Health (Academician Zhong Nanshan Research Base)


International Exhibitions/Conferences

MEDICA 2018/2019

AACC 2019

MEDLAB Middle East 2019/2020

Medicall 2019(India)

MEDLAB Asia 2019

INTRODUCTION

- The novel coronaviruses belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.
- The novel coronavirus invades human cells by the specific binding of its spike glycoprotein (ligand) to the ACE2 receptor located on human cellular membrane. In this test the ACE2 receptor has been substituted for antibody to establish a novel ligand-receptor chromatography test kit for rapid detection of the novel coronavirus.




Figure 3. Visualization of 2019-nCoV with Transmission Electron Microscopy. Negative-stained 2019-nCoV particles are shown in Panel A, and 2019-nCoV particles in the human airway epithelial cell ultrathin sections are shown in Panel B. Arrowheads indicate extracellular virus particles, arrows indicate inclusion bodies formed by virus components, and triangles indicate cilia.

SARS-CoV-2 INFECTION TIME WINDOW

Test method	PCR	IgM / IgG antibody rapid test	Antigen rapid test (latex)
Environmental requirements	High requirements for laboratory environment	Low environmental requirements, not necessarily performed in the laboratory	Low environmental requirements, not necessarily performed in the laboratory
Operational requirements	Professional testing personnel, need equipment	The operation is simple and no instrument is needed	The operation is simple and no instrument is needed
Sample type	Nasopharynx swab, high sampling requirements, pain	Blood samples, invasive sampling	Oropharyngeal saliva, sputum and feces of the posterior oropharynx
Performance	Gold standard of diagnosis;	Auxiliary diagnosis; Good sensitivity and specificity	Auxiliary diagnosis; Higher sensitivity, 100% specificity
Detection time	2-3 Hours	15 minutes	10-15 minutes
Cost	higher	Economics	Economics
Applicability	Applicable to all periods	Weeks after symptom onset	Early detection and mass screening of suspected population or animals

Features

COVID-19 ANTIGEN RAPID TEST (Latex)

Patent Application No. (USA) : T13520.PROV


JOINSTAR COVID-19 ANTIGEN RAPID TEST has complete export qualifications; non-invasive; saliva (oropharyngeal), sputum and stool can be detected, early diagnosis reassures your mind

- **Internationally innovative**, direct detection of pathogen S protein, not affected by virus mutation, high sensitivity & specificity, and can be used for early screening;
- **Convenient and non-invasive sampling**. Specimen type: oropharyngeal saliva/sputum/stool, which can be used for home self-inspection during the quarantine, and screening before resumption of work and school; Non-invasive testing is particularly suitable for continuous monitoring of children and the elderly;
- **One-step method**, easy to operate, reducing missed or false inspections caused by operator errors;
- **No equipment required**, fast detection, results are available in 10-15 minutes;
- **Storage temperature: 2~30°C**. No cold-chain transportation needed;
- **Specification: 25 tests/box, 1 test/box; Diverse cooperation modes: OEM/ODM accepted.**

Specification	P/N	Component	
25tests/box 1000tests/CTN (40*60*37CM 15/17KG)	FLCOVA200	Test cassette*25 Sample extraction tube*25	Paper cup*25 Dropper*25 Package insert*1
1test/box 400tests/CTN (40*60*37CM 6/7KG)	FLCOVA100	Test cassette*1 Sample extraction tube*1	Paper cup*1 Dropper*1 Package insert*1


Performance Characteristics

Testing Procedure

Please read the instructions carefully and allow the test device and specimens to equilibrate to temperature (15°C -30°C) prior to testing.

I. Sputum/oropharyngeal saliva (at least 200 μ L)


1. Rinse and spit with water. This is important to make sure there won't be mouth bacteria in the sputum collected.

**It is recommended to collect Sputum/oropharyngeal saliva in early-morning (first thing in the morning upon awakening, before teeth brushing, mouth rinsing, and eating breakfast)*

2. Cough deeply, Make the noise of "Kruuuu" from the throat to clear sputum/oropharyngeal saliva from deep throat.

3. Once the sputum/oropharyngeal saliva is in your mouth, release it into the container


4. Unscrew the Sample Extraction Tube

5. Transfer **200 μ L** of fresh Sputum/oropharyngeal saliva samples from container into the Sample Extraction Tube.

6. Tighten the sample extraction tube


Testing Procedure

7. Shake and mix completely

8. Take the test cassette from the packaging bag, place it on a table

9. Cut off the protrusion of the collection tube



10. Add 3 drops of the sample into the sample hole vertically

11. Read the result after 15 minutes.
If left unread for 20 minutes or more, the results are invalid, and a repeat test is recommended.

Testing Procedure

Please read the instructions carefully and allow the test device and specimens to equilibrate to temperature (15°C-30°C) prior to testing.

II. Stool Sample (at least 30mg)





1. Unscrew the Sample Extraction Tube

2. Use the sampling rod to pick up about 30mg fresh stool samples (Cover the end of the sampling rod)

3. Tighten the sample extraction tube

4. Shake and mix completely

5. Take the test cassette from the packaging bag, place it on a table

6. Cut off the protrusion of the collection tube

7. Add 3 drops of the sample into the sample hole vertically

8. Read the result after 15 minutes. If left unread for 20 minutes or more, the results are invalid, and a repeat test is recommended.

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DECLARATION OF CONFORMITY

Manufacturer: Joinstar Biomedical Technology Co.,Ltd.

Address: 10th Floor, Administration Building, NO.519, XingGuo RD., Yuhang Economic and Technological Development Zone, Hangzhou, Zhejiang, China, 311188

EC Representative's Name: Lotus NL B.V.

EC Representative's Address: Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.

Declares, that the product

Product Name and Model:
COVID-19 Antigen Rapid Test (Latex)
1 Test/Kit, 25 Tests/Kit

as described above are in conformity with the requirements as defined in IVDD98/79/EC Annex III.

Additional information:

Conformity assessment route: Directive 98/79/EC, Annex III
Classification: List Others

I, the undersigned, hereby declare that the medical devices specified above conform with the Directive 98/79/EC on in vitro diagnostic medical devices and pertinent essential requirements.

Date Signed: 2020.09.02

Xuyi ZHOU
General Manager
Joinstar Biomedical Technology Co.,Ltd.

Joinstar Biomedical Technology Co.,Ltd.

Version:0.0

Certificate

No. Q5 087635 0004 Rev. 01

Holder of Certificate: **JOINSTAR BIOMEDICAL TECHNOLOGY CO., LTD.**
10th Floor, Administration Building
No.519 Xingguo Rd.
Yuhang Economic and Technological Development Zone
311188 Hangzhou
PEOPLE'S REPUBLIC OF CHINA

Facility(ies): JOINSTAR BIOMEDICAL TECHNOLOGY CO., LTD.
10th Floor, Administration Building, No.519 Xingguo Rd., Yuhang Economic and Technological Development Zone, 311188 Hangzhou, PEOPLE'S REPUBLIC OF CHINA

Facility(ies): JOINSTAR BIOMEDICAL TECHNOLOGY CO., LTD.
No. 1 Factory Building, No. 519 Xingguo Rd., Yuhang Economic and Technological Development Zone, 311188 Hangzhou, PEOPLE'S REPUBLIC OF CHINA

Certification Mark:

Scope of Certificate: **Design, Development, Production and Distribution of Biochemical Reagent, ELISA Reagent, Clinical Laboratory Instruments and Rapid Diagnostic Reagents**

Applied Standard(s): EN ISO 13485:2016
Medical devices - Quality management systems - Requirements for regulatory purposes (ISO 13485:2016)
DIN EN ISO 13485:2016

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). See also notes overleaf.

Report No.: SH2087401

Valid from: 2020-05-27

Valid until: 2023-05-26

Date, 2020-05-07

Christoph Dicks
Head of Certification/Notified Body

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证书编号: 20200007
Certificate NO.: 20200007

产品名称: 见附件 (共1页)
Product(s): See Attachment (1 Page)

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Manufacturer: Joinstar Biomedical Technology Co., Ltd.

生产企业住所: 浙江杭州余杭经济开发区兴国路 519 号
Address of manufacturer: No.519 XingguoRD, Yuhang Economic and Technological Development Zone, 311118, Hangzhou, P.R. China

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Manufacturer: Joinstar Biomedical Technology Co., Ltd.

出口企业住所: 浙江杭州余杭经济开发区兴国路 519 号
Address of manufacturer: No.519 XingguoRD, Yuhang Economic and Technological Development Zone, 311118, Hangzhou, P.R. China

兹证明上述产品未在中国注册, 尚未进入中国市场, 该产品出口不受限制。
This is to certify that the above product(s) are not registered in China and not distributed on the Chinese market. The exportation of the product(s) is not restricted.

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Zhejiang Provincial Association For Medical Equipment Industry
(浙江省医疗器械行业协会)

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