COVID-19 Antigen Rapid Test Kit (Saliva&Swab)


Product Catalogue

COVID-19 Antigen Rapid Test Kit(Saliva/Swabs)

Funnel saliva collection kit
Selective option

Wuhan EasyDiagnosis Biomedicine Co., Ltd.

- **TOP 3** in China POCT market
 - Official Covid-19 Test Kit Supplier
- 15%** of the Market Share in China

Wuhan EasyDiagnosis Biomedicine Co., Ltd.

- Founded in **2008**
- Headquarters **Wuhan**, China
- Listed on the Shenzhen Stock Exchange in **2018**
- Employees **1000+** Globally-Based
- Export to **100+** Countries and Regions


Invest **20%** of
the main business
revenue

Patents
100+

Master or PHD
70%

New products
5 annually

Production Capacity

- Cleanroom workshop area **20,000+m²**
- **ISO13485** certificate
- Daily Production Capacity of COVID-19 Antigen Test Kit Over **2M Tests/Day**

- Daily Production Capacity of COVID-19 Nucleic Acid Test Kit Over **6M Tests/Day**

Domestic Market

- **10,000+** Hospitals, Medical Institutions and Laboratories are Using Easydiagnosis Products
- Official COVID-19 Nucleic Acid Test Supplier of **25 Provinces**, Ranking the **Top 3** in China

- Our Products Has been Distributed to **100+ Countries and Regions**


Certificate

**Quality Management System
EN ISO 13485:2016**

Registration No.: SX 2055510-1

Organization: Wuhan EasyDiagnosis Biomedicine Co., Ltd.
Room 3 & 4, 2nd Floor, Bldg 25, Phase 3.1 Wuhan Optics Valley
International Biopharmaceutical Enterprise Accelerator, No.388, Gaoxin
2nd RD, East Lake Hi-Tech Development Zone
Wuhan, 430074 Hubei P.R. China

Scope: Design and Development, Manufacture and Distribution of In-Vitro-Diagnostic Test Kits for the Diagnosis of Cancer, Amniorrhexis, Thrombotic Diseases and Hypercoagulation, Cardiac Markers, Kidney Function Testing, Immune Status, Pregnancy Testing, Prostate Function, Diabetes Testing, Neurodegenerative Diseases, Endocrine Disorders, Infection Diseases, Disease Status, Blood Gases and Genetic Testing as well as for the Monitoring of the Disease Status including Near Patient/Point of Care, Immune Quantitative Analyzer, Full-automatic Chemiluminescence Analyzer, Chemiluminescence Immunoassay Analyzer, HbA1c Analyzer, Blood Gas Analyzer, Portable ECG, Fluorescence Immunity Analyzer, Real-time PCR System.

武汉明
Wuhan Eas

The Certification Body of TÜV Rheinland LGA Products GmbH certifies that the organization has established and applies a quality management system for medical devices. Proof has been furnished that the requirements specified in the abovementioned standard are fulfilled. The quality management system is subject to yearly surveillance.

Report No.: 190129510 110
Effective date: 2021-01-27
Expiry date: 2024-01-26
Issue date: 2021-01-25

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COVID-19(SARS-CoV-2) Antigen Test Kit

Faster Detection for Earlier
Diagnosis of COVID-19

Specificity | Accuracy | Sensitivity
99.26 % | **98.50 %** | **96.15 %**

Fast

The Result Is Available
In 15 Minutes

Early

Identify Acute or
Early Infection


Convenient

No Need Instrument and
Fewer Steps Needed to
Get Result


Comparison of Different Test Technologies

	Antibody Test	Molecular Test	Antigen Test
What to Test	COVID-19 IgG/IgM Antibody	RNA of SARS-CoV-2	Protein of SARS-CoV-2
Ideal Time to Take Test	5-7 Days after Onset	3-7 Days after Infection	3-7 Days after Infection
Time to Get Result	<20 Minutes	>1.5 Hours	15 Minutes
Operation	Simple	Complex	Simple

Measurement Procedure

Interpretation of Result

Main Components

Components	Storage Temperature
Test Cassette	2°C~30°C
Antigen Extract R1	
Saliva Collector (Paper Bag or Funnel collector)	Room Temperature
Pipette /Dropper	
Extraction Tube	
Tube Stand	

0086-27-
65523649

Wuhan EasyDiagnosis Biomedicine Co., Ltd.
www.easydiagnosis.com.cn sales@easydiagnosis.com.cn

Room 3 & 4, 2nd Floor, Bldg 25, Phase 3.1, Wuhan Optics Valley International Biopharmaceutical Enterprise Accelerator,
 No. 388, Gaoxin 2nd Road, East Lake Hi-Tech Development Zone, Wuhan, 430074, Hubei, P.R. China

COVID-19 Antigen Rapid Test Kit (Saliva/Swabs)

Instruction for Use

[Product name]

COVID-19 Antigen Rapid Test Kit (Saliva/Swabs)

[Packing specification]

Saliva:

20 Tests/Kit, 10 Tests/Kit, 5 Tests/Kit, 2 Tests/Kit, 1 Test/Kit

Oropharyngeal/Nasopharyngeal swabs:

40 Tests/Kit, 20 Tests/Kit, 5 Tests/Kit, 2 Tests/Kit, 1 Test/Kit

[Intended use]

This product is intended for in vitro qualitative detection of the nucleocapsid protein antigen from SARS-CoV-2 in saliva and oropharyngeal (throat)/nasopharyngeal swabs from individuals suspected of COVID-19 by their healthcare provider. The kit is intended for use by laboratory trained personnel.

Coronaviruses are a large family of viruses which could cause illness in animals or humans. The Novel Coronaviruses (SARS-CoV-2) belong to the RNA virus of the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

[Test principle]

This kit employs immunochromatography for detection. The

specimen will move forward along the test card under capillary action. If the SARS-CoV-2 nucleocapsid protein antigen is present, they will be bound to the colloidal gold-labeled SARS-CoV-2 specific antibodies. The immune complex will be captured by coronavirus monoclonal antibody fixed in the T line. A fuchsia line would form, and the test result would be positive. If the line does not show color, the negative result will be displayed. The test card also contains a quality control line C, which shall appear fuchsia regardless of whether there is a detection line.

[Main components]

Saliva:

Spec. Components	20 tests/kit	10 tests/kit	5 tests/kit	2 tests/kit	1 test/kit
Antigen test cassette	20 pieces	10 pieces	5 pieces	2 pieces	1 piece
Desiccant	20 pieces	10 pieces	5 pieces	2 pieces	1 piece
Antigen extract R1	0.5mL X20 tubes	0.5mL X10 tubes	0.5mL X5 tubes	0.5mL X2 tubes	0.5mL X1 tube
Instruction for use	1 piece	1 piece	1 piece	1 piece	1 piece
Saliva Collector	20 pieces	10 pieces	5 pieces	2 pieces	1 piece
Saliva Dropper	20 pieces	10 pieces	5 pieces	2 pieces	1 piece
Antigen extraction tube (with tube and dropper head)	20 pieces	10 pieces	5 pieces	2 pieces	1 piece
Tube stand	1 piece	1 piece	1 piece	1 piece	1 piece

Wuhan EasyDiagnosis Biomedicine Co., Ltd.

Oropharyngeal/Nasopharyngeal swabs:

Spec. Components	40 tests/kit	20 tests/kit	5 tests/kit	2 tests/kit	1 test/kit
Antigen Test cassette	40 pieces	20 pieces	5 pieces	2 pieces	1 piece
Desiccant	40 pieces	20 pieces	5 pieces	2 pieces	1 piece
Antigen extract R1	0.5mL X40 tubes	0.5mL X20 tubes	0.5mL X5 tubes	0.5mL X2 tubes	0.5mL X1 tube
Instruction for use	1 piece	1 piece	1 piece	1 piece	1 piece
Oropharyngeal/Nasopharyngeal swabs	40 pieces	20 pieces	5 pieces	2 pieces	1 piece
Antigen extraction tube (with tube and dropper head)	40 pieces	20 pieces	5 pieces	2 pieces	1 piece
Tube stand	1 piece	1 piece	1 piece	1 piece	1 piece

[Storage conditions & period of validity]


1. Store at 2°C~30°C, and it is valid for 24 months. DO NOT FREEZE.
2. After the aluminum foil bag is unsealed, the test card should be used as soon as possible .
3. For best performance, it is recommended to test immediately after sample collection.
4. The saliva sample/Oropharyngeal (throat) swabs/ Nasopharyngeal swabs in the airtight container can be stored at room temperature for 1 hours or at 2-8°C for 24 hours or at -20°C for 7 days.

[Specimen request]

Saliva:

Relax the cheeks and gently massage for 15-30 seconds with the fingers before the specimen collection. Hold the tongue against the upper and lower jaws including the roots to enrich the saliva. To collect the saliva, install the saliva collector over the collection tube or open the collection bag, and collect saliva specimen according to the following procedure:

Bring close to the lips and gently spit in the saliva collector or collection bag and let the saliva flow into the collection tube or collection bag.


● **Oropharyngeal (throat) swabs:**

Let the patient's head tilt slightly, mouth open, and make "ah" sounds, exposing the pharyngeal tonsils on both sides. Hold the swab and wipe on both sides of the patient's pharyngeal tonsils, wipe it back and forth for at least 3 times and do it with moderate pressure applied. Place the swab specimen in the extraction tube with the Antigen extract R1 added in advance, rotate the swab for about 10 seconds, and press the swab head against the wall to

release the antigen in the swab.

● **Nasopharyngeal swabs:**

Let the patient's head relax naturally, and slowly rotate the swab against the wall of the nostril into the patient's nostril to the nasal palate, and then slowly remove it while wiping. Using the same swab, wipe the other nostril in the same way; place the swab specimen in the extraction tube with the Antigen extract R1 added in advance, rotate the swab for about 10 seconds, and then squeeze the swab head from the tube wall to release the antigen in the swab

[Test methods]

1. Unseal the package and take out the antigen test cassette.
2. Place the extraction tube on the tube stand. Place the Antigen extract R1 bottle vertically downward to allow the solution to drip freely into the extraction tube without touching the edge of the tube. Add 6 drops of R1 to the extraction tube.
3. Specimen Preparation

Wuhan EasyDiagnosis Biomedicine Co., Ltd.


Saliva: Add 2 drops of saliva into the extraction tube by saliva dropper, shake the extraction tube vigorously to mix the saliva and the extraction buffer. Squeeze the tube about 10 times to allow a thorough mixing.

Oropharyngeal/Nasopharyngeal swabs: Put the swab specimen into the extraction tube, rotate the swab for about 10 seconds, and press the swab head against the tube wall to release the antigen in the swab. Squeeze the swab over the head while taking the swab out of the extraction tube to remove as much liquid as possible from the swab. Dispose of swabs according to biohazard waste disposal method.

4. Install the dropper head on the extraction tube, put 3 drops into the specimen well of the test cassette, and start the timer.
5. Read the results within 15 minutes. Strong positive results can be reported within 15 minutes, however, negative results must be reported after 15 minutes, and the results after 25 minutes are no longer valid.

[Interpretation of test results]

- **Positive result:** if both the quality control line C and the detection line appear, novel coronavirus antigen has been detected and the result is positive for antigen.
- **Negative result:** if there is only a quality control line C, the detection line is colorless, indicating that novel coronavirus antigen has not been detected and the result is negative.
- **Invalid result:** if the quality control line C is not observed, it will be invalid regardless of whether there is detection line (as shown in the figure below), and the test shall be repeated.

4. Precision: Test the precision references, the test results shall be positive with uniform color.

5. Analytical Specificity

1) Cross-reactivity:

No false positive test results for COVID-19 (SARS-CoV-2) Antigen Test Kit were observed on specimens from the following disease states or specific conditions:

Staphylococcus aureus, streptococcus pneumonia, measles, mumps virus, Adenovirus (type 3,C1,71), Mycoplasma pneumonia, Parainfluenza virus (1-4)(2), Mycobacterium tuberculosis,CoronavirusOC43, Coronavirus229E, CoronavirusNL63, CoronavirusHKU1, bordetella pertussis, Influenza B Virus (Victoria), Influenza B Virus (Yamagata), H1N1, H3N2, EBV, Coxsackievirus A16 (CVA16), Rhinovirus, Respiratory syncytial virus, Chlamydia pneumoniae, Haemophilus influenzae, Legionella pneumophila, Mycobacterium tuberculosis, Streptococcus pyrogenes, Pneumocystis jirovecii (PJP), Pooled human nasal wash .

2) Interference:

No interference was observed with the potentially interfering substances listed below at the indicated concentration:

Commonly used drugs, i.e., 5mg/mL Phenylephrine, 5mg/mL Oxymetazoline, 0.9% sodium chloride, 200µg/L beclomethasone, 150µg/L dexamethasone, 400µg/L flunisolide, 60mg/mL triamcinolone acetonide, 2mg/mL budesonide, 5mg/mL mometasone, 250µg/L fluticasone, 1mg/mL Histamine hydrochloride, 60µg/L alpha- interferon, 142ng/mL zanamivir, 0.25mg/L ribavirin, 24.6mg/mL oseltamivir, 0.1mg/L peramivir, 4mg/L Lopinavir, 250mg/L Ritonavir, 100mg/L Arbidol, 500mg/L levofloxacin, 500mg/L azithromycin, 100mg/mL ceftriaxone, 50mg/mL meropenem, 4µg/mL tobramycin, 0.9mg/ml mucin, 50% blood

Wuhan EasyDiagnosis Biomedicine Co., Ltd.

(human), Human Anti-mouse Antibody (HAMA), biotin have no effect on the test results of this kit.

6. Clinical Performance

1) Oropharyngeal swabs and nasopharyngeal swabs specimens:

Contrast Results Statistics of Clinically Confirmed/Excluded Results

Evaluation Reagent	Clinical Confirmed/Excluded Results		Total
	Confirmed	Excluded	
Positive	125	3	128
Negative	5	401	406
Total	130	404	534

Result calculation:

(1) Clinical sensitivity: 96.15%, 95% confidence interval: [91.31%, 98.35%].

(2) Clinical specificity: 99.26%, 95% confidence interval: [97.84%, 99.75%].

(3) Clinical accuracy: 98.50%, 95% confidence interval: [97.07%, 99.24%].

*In the stratified statistics of different stages of the disease, 52 specimens from 0-3 days, and the positive detection rate is 98.08%.

2) Saliva specimens:

Contrast Results Statistics of Clinically Confirmed/Excluded Results

[Limitations of inspection methods]

1. This reagent is for in vitro diagnostic use only.
2. This reagent is only used to detect human sterile swab extracts or saliva sample. The results of other specimens maybe inaccurate.
3. This reagent is only used for qualitative detection and can not indicate the level of novel coronavirus antigen in the specimen.
4. This reagent is only a clinical auxiliary diagnostic tool. If the result is positive, it is recommended to use other methods for further examination in time and the doctor’s diagnosis shall prevail.

[Product Performance Indicators]

1. **LoD:** The LoD for direct swab was established using heat-inactivated SARS-CoV-2. The strain was spiked into the pooled human oropharyngeal swab matrix obtained from multiple healthy volunteers eluted in R1 and confirmed as SARS-CoV-2negative by RT-PCR to prepare positive samples. The estimated LoD found from the initial serial dilution test was confirmed by testing 20 replicates. The confirmed LoD for direct swab was 5 x 10² TCID₅₀/ml.
2. **Negative Agreement:** Test the negative reference, and the negative accordance rate shall be 100%.
3. **Positive Agreement:** Test the positive reference, and the positive accordance rate shall be 100%.

Evaluation Reagent	Clinical Confirmed/Excluded Results		Total
	Confirmed	Excluded	
Positive	123	4	127
Negative	5	423	428
Total	128	427	555

Result calculation:

(1) Clinical sensitivity: 96.09%, 95% confidence interval: [91.49%,98.01%].

(2) Clinical specificity: 99.06%, 95% confidence interval: [97.72%, 99.53%].

(3) Clinical accuracy: 98.38%, 95% confidence interval: [97.00%, 99.09%].


*In the stratified statistics of different stages of the disease, 52 specimens from 0-3 days, and the positive detection rate is 98.04%.

[Precautions]

1. This reagent must be used by trained or professional clinical testing personnel by following all laboratory management regulations.
2. Read the instructions carefully before using the kit, and strictly control the reaction time. If you do not follow the instructions, you will get inaccurate results.
3. The specimen shall be tested in a laboratory with certain conditions. All specimens and materials during testing should be handled in accordance with the laboratory practice for infectious diseases.
4. Keep away from moisture, do not open the aluminum platinum bag before it is ready for testing. Do not use when the aluminum foil bag is damaged or the test card is damp.
5. Please use it within the validity period.
6. Wait all reagents and specimens come to room temperature (15 ~ 30 °C) before use.

7. Do not replace the components in this kit with components in other kits.
8. Do not dilute the specimen for testing, otherwise you may get inaccurate results.
9. The kit shall be stored in strict accordance with the conditions specified in this manual. Please do not store the kit under freezing conditions.
10. The test methods and results must be interpreted in strict accordance with this specification.

[Index of Symbols]

	Temperature Limit		Use-by date
	Batch/Lot code		In vitro diagnostic medical device
	Manufacturer		Catalogue number
	Contains sufficient for <N> tests		Consult instructions for use
	Do not re-use		Authorized representative in the European Community
	Date of manufacture		

Wuhan EasyDiagnosis Biomedicine Co., Ltd.

Do not use if package is damaged

Sterilized using irradiation

CE Certification

[INFORMATIONENINQUIRIES AND GENERAL INFORMATION]

Wuhan EasyDiagnosis Biomedicine Co., Ltd.

Address: Room 3 & 4, 2nd Floor, Bldg 25, Phase 3.1 Wuhan Optics Valley International Biopharmaceutical Enterprise Accelerator, No.388, Gaoxin 2nd RD, East Lake Hi-Tech Development Zone, 430074 Wuhan, China

Tel: +86(0)27-87808955

Fax: +86(0)27-87808005

WEB: www.mdeasydiagnosis.com

Email: info@ediagnosis.cn

MedNet EC-REP GmbH

Borkstrasse 10, 48163 Münster, Deutschland

[SWAB INFORMATION]

ShenzhenKangDaAn Biological Technologyco, LTD

East- 1, 3rd floor, Building 2, Shunheda factory, Liuxiandong industrial zone, Xili street, Nanshan district, Shenzhen .

Tel: +86 0755- 82836262

Fax: +86 0755- 83426595

WEB: www.kdasw.com

Name: Share Info Consultant Service LLC Repräsentanzbüro
Address: Heerdter Lohweg 83, 40549 Düsseldorf

**EG-KONFORMITÄTSERKLÄRUNG · EC DECLARATION OF CONFORMITY
DÉCLARATION CE DE CONFORMITÉ · DICHIARAZIONE CE DI CONFORMITÀ**

Name und Adresse des Herstellers: /
Name and address of the manufacturer: /
Nom et adresse du fabricant: /
Nome e indirizzo del fabbricante:

Wuhan EasyDiagnosis Biomedicine Co., Ltd
Room 3 & 4, 2nd Floor, Bldg 25, Phase 3.1 Wuhan Optics Valley
International Biopharmaceutical Enterprise Accelerator, No.388,
Gaoxin 2nd RD, East Lake Hi-Tech Development Zone,430074
Wuhan, China

EU-Vertreter/
Authorized EU Representative/
représentants européens /
Presentanti dell'UE:
DIMDI No.:

**MedNet EC-REP GmbH Borkstrasse 10, 48163 Münster,
Germany**
DE/0000048589

Wir erklären in alleiniger Verantwortung, dass / We ,as manufacturer,declare under our sole responsibility that /
Nous déclarons sous notre propre responsabilité que / Dichiariamo sotto la sola responsabilità che

das Medizinprodukt: /
the medical device: /
Product Name:/
Type/Model:
le dispositif médical: /
il dispositivo medico:

COVID-19 Antigen Rapid Test Kit (Saliva/Swabs)
**Analyte:antigen from SARS-CoV-2 in saliva and oropharyngeal
(throat)/nasopharyngeal swabs**

Saliva:
20 Tests/Kit, 10 Tests/Kit, 5 Tests/Kit, 2 Tests/Kit, 1 Test/Kit
Oropharyngeal/Nasopharyngeal swabs:
40 Tests/Kit, 20 Tests/Kit, 5 Tests/Kit, 2 Tests/Kit, 1 Test/Kit

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Others

Nach Richtlinie 98/79/EG / selon directive 98/79/CE
secondo direttiva 98/79/CE / according to direct. 98/79/EC

den einschlägigen Bestimmungen der Medizinprodukte-Richtlinie 98/79/EG und deren Umsetzungen in nationale
Gesetze entspricht.
meets the provisions of the directive 98/79/EC and its transpositions in national laws which apply to it.
remplit toutes les exigences de la directive sur les selon directive 98/79/CE et de ses transpositions en droit national
qui le concernent.
soddisfa tutte le disposizioni della direttiva 98/79/CE e della loro trasposizione nel diritto nazionale che lo riguardano.

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in its function of the European Authorized Representative, in accordance with the In Vitro Diagnostic Directive 98/79/EC, hereby confirms the submission of the registration of the following *in vitro* diagnostic medical devices into the German DIMDI data base

COVID-19 Antigen Rapid Test Kit (Saliva/Swabs)

Analyte: antigen from SARS-CoV-2 in saliva and oropharyngeal (throat)/nasopharyngeal swabs

on behalf of

Wuhan EasyDiagnosis Biomedicine Co., Ltd.

Room 3 & 4, 2nd Floor, Bldg. 25, Phase 3.1 Wuhan Optics Valley International Biopharmaceutical Enterprise Accelerator, No. 388, Gaoxin 2nd RD, East Lake Hi-Tech Development Zone, 430074 Wuhan, P.R. China

according to the directive 98/79/EC of the European Parliament and of the Council of the European Union relating to *in vitro* diagnostic medical devices.

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Formblatt für In-vitro-Diagnostika / Form for In Vitro Diagnostic Medical Devices

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Only for registration use or business use

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Bezeichnung / Name Wuhan EasyDiagnosis Biomedicine Co., Ltd.